

Et, FiAgents & Waveform of Halothane, Enflurane, Isoflurane, Sevoflurane or Desflurane

Note: Need to manually select one of the above five anesthetic gases

## Features of this Multigas Monitor

- An economic monitor for anesthetic agent, which supports safe agent administration
- Reliable and accurate measurement, based on Kingst's strong heritage technology of infrared gas analysis
- Special designed water filter can be used keeping work under the heavy humid environment, also suitable for an intubated or non-inbutared patients.
  - Audio-optical alarm for exceeding limits and Apnea
  - 24 hours trend diagram keeps still available after power-off
  - Data output via USB
  - Large Capacity lithium battery for 10 hours UPS
  - Nellcor or our Digital SpO2 and adult, pediatric or neonate.
  - The Analog Output Module is optional

## Technical Specifications

Technique: NDIR (Non-dispersive infrared gas analysis)

Power: AC100~240V, 50/60Hz; DC3.7V;

Temperature: 5 to 50° C

Humidity:<93%

Atmospheric pressure: 70-120kP

Battery: Lithium battery, 3.6V, 3800mAh, 6 hours (approximately)

Charging time: 4 hours from empty to full **Dimension:** 72mm(W)x155mm(H) x38mm(D)

Weight: 360 grams
Display: 3.5' 'TFT

CO2 Range: 0-19.7%, 0-150mmHg, or 0-20kPa
Accuracy: 0 - 40 mmHg ± 2 mmHg
41 - 70 mmHg ± 5% of reading
71 - 100 mmHg ± 8% of reading
101 - 150 mmHg ± 10% of reading

EtCO2 update: every breath

Trends: 24 hours
N2O Range: 0-100%

Accuracy: ± (2 vol% + 2% of reading)

HAL, ENF, ISO Range: 0-12%

Accuracy: 0-8%:  $\pm (0.15 \text{ vol}\% + 5\% \text{ of reading})$ 

8-12%: unspecified

SEV Range: 0-15%

Accuracy: 0-10%:  $\pm (0.15 \text{ vol}\% + 5\% \text{ of reading})$ 

10-15%:unspecified

**DES** Range: 0-25%

Accuracy: 0-18%:  $\pm (0.15 \text{ vol}\% + 5\% \text{ of reading})$ 

18-25%: unspecified

All FiAG update: every breath

All FiAG Trends: 24hours

Respiration rate Range: 3-150 bpm

Accuracy: the bigger value of 1% of reading and ± 1bpm

Trends: 24hours

Water filter Detection: Insertion automatically turns sampling pump

on. Removal automatically turns sampling pump off

Flow set: 50-150ml/min selectable

Exhaust Port: Yes

Auto-compensation for atmospheric pressure variation

Oximeter Data (optional)

Oximetry Saturation Range: 0-100%

Accuracy: ±3% when 50-69%

±2% when 70-100%

Trend length: 24hours
Pulse Rate Range: 30-250 bpm

Accuracy: the bigger value of  $\pm 2$ bpm or  $\pm 2\%$  (reading)

Trend length: 24hours

Audible and visual alarm system

Adjustable High and Low EtCO2,RR and FiAG,No Adapter,

Clogged Adapter, No Breath (Apnea), Low Battery, so on.

Regulatory: Designed to meet IEC 60601-1:2005+A1:2012(E):Medical electrical equipment – Part 1: General requirements for basic safety and essential performance; IEC 60601-1-2:2014: Medical electrical equipment – Part 1-2: General requirements for basic safety and essential performance – Collateral Standard: Electromagnetic disturbances – Requirements and tests; Broad-band random vibration test in accordance with IEC 60068-2-64:2008; Shock test in accordance with IEC 60068-2-27:2008; Free-fall to IEC 60068-2-31:2008

Designed to comply with 93/42/EEC (MDD CE Marking), FDA Standards, Minimum Performance and Safety Requirements for Capnometers ISO 80601-2-55-2011: Medical Electrical Equipment performance requirements for the basic safety and essential performance of respiratory gas monitors

No futher notification if some technical specification change. Please take update technical specification or test result as the stardand



## Accessories: Water filter, Sample lines and Airway Adap





T4F Filter 664522





Sample line 662503



Adult Nasal tube Infant Nasal tube 660424



CO2/O2 Nasal tube 660425 CO2/O2 infant Nasal tube 660426







Beijing Kingst Commercnd Trade Co. Ltd

Add: 5-6F, Building 3, No.27, Yongwang Road, Z-park Daxing Bio-medicine Industry Park, Dist. Daxing, Beijing P.R.C Code: 102600